





















SSRL Radioactive Material Sample Holder Catalog

7/25/08

Hazard Class Category Containment #	3 Layer containment for Very High Radiotoxicity (Group 1)	
1.a	<p>LBNL Lexan or aluminum sample holder with kapton tape surrounded by 2 each individual heat sealed plastic bag.</p> <p>Layer 1- Kapton Tape, sealed Layer 2- Heat sealed plastic bag Layer 3- Heat sealed plastic bag</p> <p>Physical Approvals: Ambient temperature</p>	  <p>Material: PCTFE or Lexan (polycarbonate), 1/16 in.</p>
1.b	<p>Sample holder with lexan windows and indium seam</p> <p>Layer 1-Lexan window with indium seal</p> <p>Layer 2-Aluminum with Kapton tape, screws</p> <p>Layer 3-durable plastic bag or BL11-2 prep room tent inside prep room. Cryostat S Steel housing assembly with kapton windows and bolting ring outside of BL 11-2 prep room</p> <p>Physical Approvals: Ambient temp Cold</p>	
1.c	<p>LANL cryostat sample holder</p> <p>Layer 1- Layer 2- Layer 3-</p> <p>Physical Approvals: Ambient temp Cold</p>	

<p>1.d</p>	<p>USGS cryostat holder Layer 1- Layer 2- Layer 3-</p> <p>Physical Approvals:</p>	
<p>1.e</p>	<p>G-XAS cell One sample per cell = ^{231}Pa sorbed onto TiO_2 single crystal. Maximum activity : 15000Bq - Doubly contained crystal in two sealed polyethylene envelopes. Check for no contamination of each envelope. - External envelope glued onto the cell surface. - Cell cap (polyamide 0.7 mm thick) sealed with Viton o-ring. - Additional glue in the thread. Layer 1-polyethylene Layer 2- Layer 3- Cell cap (polyamide 0.7 mm thick) sealed with Viton o-ring.</p> <p>Physical Approvals: Nominal operating conditions of BL 11-2 at 17 keV (focused beam).</p>	<p>G-XAS cell for radioactive samples, BL-11-2 $^{231}\text{Pa}/\text{TiO}_2$, 15 000 Bq</p> 
<p>1.f</p>	<p>Diamond anvil cell Layer 1 Layer 2- Layer 3-</p> <p>Physical Approval: Pressure</p>	
<p>1g</p>	<p>LBNL Triple contained aluminum holders with Kapton or Mylar windows. Fit into cryostats. Layer 1-nested, 1 side is epoxied window, other side is indium wire pressed with lid and cap screws. Layer 2 – nested, 1 side is epoxied window, other side is indium wire pressed with lid and cap screws. Layer 3-Outer 1 side is epoxied window, other side is indium wire pressed with lid and cap</p>	

	<p>screws.</p> <p>Physical Approvals: Ambient temp, Cold</p>	
1h.	<p>LANL Conradson</p> <p>Layer 1 polystyrene Layer 2 Stainless Steel Chamber Layer 3 Glovebag</p>	  
	<p>1 Layer containment for Low Radiotoxicity (Group 4)</p>	
4.a		
4.b	<p>Slotted metal sample holder with unpolished polystyrene sample material Layer 1-polystyrene</p>	 
4.c	<p>Slotted metal Al sample holder with kapton tape and bolting ring Layer 1-Epoxyed sealed kapton tape 1 side, indium sealed cap screw inside bolting ring, cap screws on kapton tape other side.</p> <p>Physical approvals: Room temperature</p>	
4.d	<p>Epi tube, screw top Layer 1-plastic vial, screw top with rubber gasket. Plastic seal bag secondary.</p> <p>Uses-Liquid sample in 2 mL, screw-top, polypropylene centrifuge tube inside three layers of heat-sealed plastic.</p>	 

	<p>Physical approvals: Room temperature</p>	
<p>4.e</p>	<p>Epi tube, heat sealed Layer 1-Sample tube with glue? sealed cap in plastic bag.</p> <p>Uses: Lukins, LBNL Powder sample in heat-sealed plastic tube (actually a plastic pipette in this case although I also use heat sealed Epi tubes) inside three layers of heat-sealed plastic.</p> <p>Physical approvals: Room temperature</p>	
<p>4.f</p>	<p>Layer 1-quartz capillaries with flame sealed ends. X-ray scattering. Layer 2-Must be used with a secondary protective containment.</p> <p>Physical approvals: Room temperature</p>	
<p>4.g</p>	<p>John Barger Holder kapton (polyimide) or mylar (PET, polyester) or mica window adhesive-backed film with thickness .0025 to .010 in. Windows will be self-secured (via adhesive backing) to sample holders. A layer of window material will be applied to overlap around the margins of the sample holder to reinforce the seal provided by the adhesive of the primary tape layer.</p> <p>Layer 1 –Capton , Myl,ar Mica window Al gasket ID of bolting ring</p>	

<p>4.h</p>	<p>Bolted Aluminum Holder Layer 1- Kapton Tape Aluminum gasket</p> <p>Physical approvals: - Room temperature - Cold</p>	
<p>4.i</p>	<p>Carbon Films, LANL loaded inside vacuum chamber device Uranium Layer 1- Polished polystyrene inside plastic. (polished polystyrene alone is not approved as primary containment)</p> <p>Physical Approvals:</p>	
<p>4j</p>	<p>LBL Kapton film with indium seal and cap screws.</p> <p>Physical approvals Ambient, Cold, vacuum</p>	
<p>4k</p>	<p>SLAC, Uranium capillary</p> <p>Kapton film with plexiglass and aluminum/plastic caps. capscrews and O- ring on ends.</p>	