

A Quick Guide to XAS

This quick-guide is intended as a general overview of procedures for successfully performing X-ray Absorption Spectroscopy experiments.

In general, a good XAS experiment will be guided by the **HALO** principle:

- **Harmonics**
 - reduce energy harmonics by properly tuning monochromator and/or using harmonic rejection mirror
- **Alignment**
 - align mono, slits, detectors, and sample
- **Linearity**
 - operate ion chambers, Lytle detector, and Germanium detector in a linear response regime.
- **Offsets**
 - take offsets after any change to detector electronics and before collecting data

XAS Support Staff

XAS support cell phone

- Mon – Fri: 7 am – 10 pm
- Sat – Sun: 8 am – 5 pm
- Matthew Latimer 510-259-8102
- Erik Nelson 415-218-5505

beam-line optimization, detectors, cryostats, data collection

- Matthew Latimer X4944
- Erik Nelson X5103

equipment: rails, cryostats, liquid helium

- Jeff Maske x 3833

Ionization Chambers

high voltage

- check that detector is **not** space charge limited
 - check at beginning-of-scan energy, fully tuned
 - if space charge limited, may be tricky to tell if fully tuned
 - increasing high voltage should **not** increase I0
 - check when changing energy regime, especially when changing gas

gas

- select gas based on energy
 - energy < 5 keV ⇒ helium
 - 5 keV < energy < 15 keV ⇒ nitrogen
 - energy > 15 keV ⇒ argon
- set flow to a few bubbles per second (do **not** leave bubbler bubbling)

amplifiers

- if using SR570 amps
 - run “XASUTILS” program and select “ION CHAMBER CONTROLS”
 - under Controls / More Controls, check filter settings
 - filter type = low pass 6 dB
 - low frequency = 30 Hz
 - set gain such that beam-on gives ≈ 1 (< 5)
 - check value at beginning-of-scan energy
 - check gain after changing offset
 - set offset such that beam-off gives ≈ 0.004 (> 0)
 - check offset after changing gain
- if using Keithley preamps
 - set Rise Time on amplifiers to 30 (Lytle) or 10 (Germanium), 10 is a good all-purpose setting
 - set Gain such that beam-on gives 1 to 5 (< 10)
 - check value at beginning-of-scan energy
 - check gain after changing offset
 - set Suppression on inverse setting of Gain
 - adjust Zeroing Knob such that beam-off gives ≈ 0.004 (> 0)
 - check offset after changing gain

Lytle Detector

gas

- select gas based on fluorescence energy
 - fluorescence energy < 3 keV ⇒ nitrogen
 - 3 keV < fluorescence energy < 7 keV ⇒ argon
 - fluorescence energy > 7 keV ⇒ krypton (ideal), argon (< 12 keV)
- gas change procedure: **use care** – windows are fragile
 - 1) start with gas tubing disconnected from detector and gas flow set to 0
 - 2) slowly turn up gas flow to a barely detectable level
 - 3) connect gas tubing to detector: source to bottom port, return to top port

Soller slits and filter

- can use filter to reduce scatter signal
 - in general, for measuring $K\alpha$ radiation of element with atomic number Z , use filter with atomic number $Z-1$
 - $K\beta$ of element $Z-1$ will overlap $K\alpha$ of element Z
 - is the sample's fluorescence signal greater than the filter's?
- Soller slits reduce fluorescence from filter

preamplifier

- set gain such that beam on gives 1 to 5 (< 10), but avoid gain of 1000
 - check value above fluorescence-edge-energy
 - check gain after changing offset
- set offset such that beam off gives ≈ 0.004 (> 0)
 - adjust using potentiometer R9
 - check offset after changing gain
- Rise Time is 45 ms, so set Keithleys to < 30 ms or SR570s to > 10 Hz

alignment

- 1) set MONO above fluorescence-edge-energy
- 2) scan Lytle position and maximize Lytle signal
 - on 2-3, 4-1, 4-3
 - a) visually align horizontally
 - b) scan VERT
 - on 7-3 and 9-3
 - a) scan VERT
 - b) scan HOR
 - c) scan VERT again

Germanium Detector

moving detector

- **Be Careful!** - be inside the hutch and prepared to abort the move if necessary

count rates

- operate in linear response regime if not dead-time correcting
 - maximum Incoming Count Rate (ICR) versus shaping time
 - 0.125 μ sec \Rightarrow 110,000 counts/sec
 - 0.250 μ sec \Rightarrow 70,000 counts/sec
 - 0.500 μ sec \Rightarrow 40,000 counts/sec
 - maximum Single Channel Analyzer (SCA) count rate versus shaping time
 - 0.125 μ sec \Rightarrow 30,000 counts/sec
- 30 element detector anomalous behavior
 - ICRs > 200,000 counts/sec may cause gain-shifting
 - if observed, please contact XAS staff

windowing

- 1) run "GE13" or "GE30" program, depending on detector
- 2) insert calibration sample into beam
- 3) check counts
- 4) set SCA window for one detector element
 - a) activate button for selected SCA window
 - b) select Anticoincidence Gating Mode
 - c) set SCA width and centroid position in energy
 - d) adjust SCA window position: optimal width and centroid depends on effect of scatter on pre-edge
- 5) set SCA windows for remaining detector elements
 - a) after adjusting first window, select APPLY TO ALL
 - b) adjust window for each detector element, but do **not** select APPLY TO ALL afterwards


Soller slits and filter

- can use filter to reduce scatter signal
 - in general, for measuring $K\alpha$ radiation of element with atomic number Z, use filter with atomic number Z-1
 - $K\beta$ of element Z-1 will overlap $K\alpha$ of element Z
 - is the sample's fluorescence signal greater than the filter's?
- Soller slits reduce fluorescence from filter
- align using Lytle detector

Germanium Detector (continued)

Beryllium window

- removable outer window protects the thin inner window
 - outer, removable Be window is 0.25 mm
 - inner Be window:
 - Systems 1, 2, 3: 0.05 mm (0.002")
 - Systems 4, 5: 0.025 mm (0.001")
- use the plot below to decide when to remove the window (for example, at Mn ~11% of the fluorescence signal is absorbed by the outer window)
- please contact XAS staff to remove window



Cryostats

Please refer to “User Guide: SSRL Biotechnology Continuous Flow Liquid Helium Cryostats” manual for background and details.

cool-down

- 1) start with the sample-port under vacuum
- 2) close sample-port valve on cryostat (otherwise, clog may form)
- 3) set Heater mode on temperature controller to MANual
- 4) set Heater power on temperature controller to 0
- 5) increase Gas Flow Pump's Exhaust to ≈ 2 liter/hour
- 6) as desired temperature is approached, stabilize the temperature
 - adjust the Gas Flow Pump's Exhaust
 - may need to set Heater mode on temperature controller to AUTO and set SET temperature

sample change

- 1) open helium gas toggle valve
- 2) open sample-port valve on cryostat
- 3) check for gas flow from relief valve on cryostat
- 4) change samples
- 5) cycle between vacuum and helium gas five times, ending with helium
- 6) check for gas flow from relief valve on cryostat
- 7) close sample-port valve on cryostat
- 8) close helium gas toggle valve

warm-up

- 1) cap Germanium detector and retract
- 2) open vacuum toggle valve
- 3) open sample-port valve on cryostat (otherwise, explosion may result)
- 4) decrease Gas Flow Pump's Exhaust to ≈ 0.2 liter/hour
- 5) set Heater mode on temperature controller to MANual
- 6) for quick warm-up, set Heater power on temperature controller to 15
- 7) as desired temperature is approached, stabilize the temperature
 - set Heater power on temperature controller to 0
 - adjust the Gas Flow Pump's Exhaust
 - may need to set Heater mode on temperature controller to AUTO and set SET temperature

Cryostats (continued)

stuck sample rod

- 1) do **not** force the sample rod
- 2) cap Germanium detector and retract
- 3) warm-up to 45 to 65 K (see warm-up procedure above)
- 4) **slow** warm-up above 85 K (reduce Heater power in temperature controller)
- 5) **gently** remove sample rod

alignment

- 1) scan CRYOVERT
 - protein cell \Rightarrow minimize I1
 - powder cell \Rightarrow maximize I1
- 2) adjust CRYOHOR
 - fluorescence screen with camera is easier
 - film or “burn paper” is more precise

temperature range considerations

- photoreduction
- Debye-Waller effect
- sample’s state and integrity

liquid helium usage

- keep a log
 - 1) record “gas flow pump exhaust” rate (liter / hour)
 - 2) calculate usage: (time) x (exhaust rate)

Data Collection

offsets

- take offsets after any change to detector electronics
 - do before collecting data
 - under Detectors : Current, select Offsets, collect for 10 seconds

crystal glitches

- glitch library at <http://www-ssrl.slac.stanford.edu/smbin/dataextractnew.pl> , Beamlines, Interactive Crystal Glitch Library
- determine glitch's effect on EXAFS
 - run it: check for spike in data multiplied by k^3
 - if necessary, shorten EXAFS region

signal-to-noise (“Am I done with this sample?”)

- always at least two sweeps: do they agree?
- visually check on signal-to-noise of sweep
- more quantitative check:
 - 1) start with sum of sweeps
 - 2) do background subtraction
 - 3) multiply by k^3
 - 4) calculate Fourier transform
 - 5) compare successive averages and make judgment on improvement with number of sweeps
 - 6) estimate signal and noise amplitudes

Checks

daily

- liquid helium dewar level for cryostat
- liquid nitrogen level for Germanium detector dewar
- gas cylinders for the ionization chambers and Lytle detector
- He gas for cryostat sample space

after each beam fill

- mono slit vertical translation
- table vertical translation

after each mono-region change

- mono calibration
- table alignment
- gains and offsets

after each sample change

- if using Germanium detector - count rates

after each sweep

- monochromator detuning (the higher the energy, the more necessary)

always!

- quality of the data - look at it!
- cryostat temperature