

# Surfaces, interfaces, thin films & single crystal diffraction

Mike Toney (SSRL),

many others: Bao, Bargar, Brown, Chang, Chabynec, Chianelli, Clemens, Fisher, Fuller, Frisbie, Gilbert, Greven, Lucas, Leisch, McGehee, McIntyre, Nilsson, Pelton, Salleo, Subramanian, Suzuki, Waychunas



1. Surface & Interface, Thin Film & Single Crystal Science
2. Beam Line Needs and a Solution

# Surface & Interface Science



Surfaces and interfaces widespread importance:

- catalysis, e.g., fuel cells (Toney, Leisch, Nilsson, Chianelli)
- minerals and oxides (Brown, Bargar, Waychunas, Gilbert)
- electrochemistry (Lucas, Toney, Leisch)



What is interface atomic structure and how does this affect reactivity?

Need:




- intensity (brightness)
- flexible manipulation of sample orientation
- ability to hold cells and sample chambers
- rapid data acquisition

# Thin Film Science



Thin films widespread importance:

- organics, e.g, for transistors, solar cells (Chabinye, Bao, McGehee, Salleo, Subramanian, Frisbie, Fuller)
- inorganics, e.g, electronics, magnetics (Chang, McIntyre, Clemens, Suzuki, Salleo, Pelton)



What is atomic/molecular structure of film and how does this affect properties?

Need:



- intensity (brightness)
- ability to tune penetration depth
- small  $\alpha$  but large  $\beta$  (Qz)
- simultaneous Qz data collection
- ability to hold environmental chambers

# Single Crystal Science



Single crystal properties important:

- physics, e.g. high  $T_c$ , phase transitions (Fisher, Greven)



Need:

- intensity (brightness)
- high resolution
- ability to hold cryostats & furnaces without interference from diffractometer

Understand physical  
(electronic, superconducting)  
properties

# High Resolution Diffractometer



- high brightness
- approximately 5-30 keV
- BL 7-2
- Psi/8-circle
  - accurate sample alignment
  - flexible detector scans
  - ability to hold dispex & other reactors & chambers
  - different detectors:
    - crystal analyzer
    - fast 1D PSD
    - small area detector